

ELEVA® Fine

Hyaluronic acid dermal filler

ELEVA® Fine is engineered as the premier solution for fine corrections, setting new standards in the management of superficial lines. This filler is formulated with a precise concentration of 22.0 mg/mL of cross-linked hyaluronic acid and is infused with 0.3% lidocaine hydrochloride, ensuring a treatment experience characterized by minimal discomfort and maximum patient comfort. The hallmark of ELEVA® Fine's engineering lies in its unique rheological profile. The combination of a low G' (modulus of elasticity) and superior spreadability allows the product to integrate seamlessly into sensitive and superficial tissues. This ensures a smooth, uniform, and completely natural aesthetic outcome, effectively mitigating the risk of lumps, bumps, or asymmetry. Each unit is supplied in a 1.1 mL ready-to-use syringe, pre-filled with a physiological buffer solution, and comes complete with two specialized 30G needles. This thoughtful design empowers aesthetic professionals with the high degree of control necessary to perform delicate corrections with unparalleled precision, even in the most sensitive facial areas.

Product information

specifications

The specialized structural design of ELEVA® Fine makes it the ideal solution for delicate facial corrections, where the convergence of science and artistry is essential for natural results. By specifically targeting superficial and sensitive skin layers, ELEVA® Fine delivers consistent, high-precision clinical outcomes in the following areas:

- Periorbital Region: Effective softening of crow's feet and specialized correction of the tear trough. Thanks to its controlled rheological profile, ELEVA® Fine prevents unwanted water absorption, delivering a refreshed and smooth appearance without the risk of puffiness or edema.

- Perioral Region: Precise correction of fine vertical lines (smoker's lines), creating a uniform surface that integrates seamlessly with the natural skin texture.
- Forehead and Surface Texture: Smoothing of fine skin wrinkles and enhancement of overall skin surface quality, resulting in a clearer, smoother, and more youthful complexion.

In essence, ELEVA® Fine is more than a dermal filler; it is a precision instrument for restoring facial harmony. By ensuring natural-looking results without the risk of unwanted irregularities, it serves as an unrivaled choice for aesthetic professionals who prioritize delicate refinement over aggressive volumization.

contraindications

ELEVA® must NOT be used in the following conditions:

- Intravascular injection (risk of vascular occlusion, embolization, ischemia, or necrosis).
- Glabellar and periorbital/tear trough injection with ELEVA® Shape and ELEVA® Contour (reserved exclusively for Fine/Balance where clinically indicated).
- Breast augmentation, gluteal enlargement, or implantation into bones, muscles, ligaments, or joints.
- Patients with known hypersensitivity to Hyaluronic Acid, Gram-positive bacterial proteins (Streptococcus), or Lidocaine/amide-type local anesthetics.
- Patients with active cutaneous inflammation or infection (acne, herpes simplex, dermatitis) in or near the target site.
- Patients with history of autoimmune disease, severe allergies, anaphylaxis, porphyria, or untreated epilepsy.
- Patients prone to keloids, hypertrophic scarring, or severe pigmentation disorders.
- Patients with coagulation disorders or receiving active anticoagulant/thrombolytic therapy.
- Pregnancy and breastfeeding.
- Patients under 18 years of age.
- Concomitant use with laser therapy, deep chemical peeling, or dermabrasion procedures until full re-epithelialization is achieved.

PRECAUTIONS & DRUG INTERACTIONS

- Medical Training: Must be administered exclusively by licensed, authorized, and certified medical practitioners trained in advanced injection anatomy.
- Incompatibilities & Chemical Interactions: Hyaluronic Acid precipitates in the presence of quaternary ammonium salts (such as benzalkonium chloride or cetylpyridinium chloride). Avoid contact between ELEVA® and instruments cleaned with the mentioned solutions.
- Drug Interactions: Patients taking antiplatelet agents, NSAIDs, high-dose Vitamin E, Vitamin C, Omega-3 fatty acids, or herbal supplements (Ginkgo Biloba, St. John's Wort) may experience increased bruising and bleeding. Caution is advised in patients taking medications that inhibit hepatic metabolism or in those with pre-existing cardiac conduction disorders (due to lidocaine content).
- Pre-existing Implants: Do not inject into anatomical sites with permanent (silicone, PMMA) or semi-permanent (e.g., PLLA) fillers.
- Magnetic Resonance Imaging (MRI): ELEVA® is MRI-safe and non-paramagnetic.

- **Product Inspection:** Inspect syringe integrity prior to use. Do not use if the sterile barrier is compromised, the gel contains turbidity or particulate matter, or if the plunger rod does not glide smoothly. Do not re-sterilize or reuse.

PRE-INJECTION PATIENT CARE PROTOCOL

- **1 Week Prior:** Discontinue elective blood thinners, NSAIDs, aspirin, high-dose Vitamin E/C, Ginkgo Biloba, and alcohol under physician supervision. Avoid tanning beds and intensive UV exposure.
- **3 Weeks Prior/After:** Avoid invasive dental procedures, oral surgery, and facial tattooing.
- **Herpes Prophylaxis:** In patients with a history of recurrent herpes labialis undergoing perioral/lip injections, prophylactic oral acyclovir is strongly recommended.
- **Skin Preparation:** Thoroughly cleanse and disinfect the skin using chlorhexidine or alcohol-based antiseptics (avoiding quaternary ammonium compounds). Topical anesthetic cream may be applied 30 minutes prior if desired.

Mavara Salamat Yas

Unit 4, No. 2, Sharifi Alley, Sanjabi St., Madar Sq., Mirdamad Blvd., Tehran, Iran

021-91301406

info@msypharmed.com

msypharmed.com



Full, up-to-date product page